

# HF100+ ultra light

# MINXRAY<sup>®</sup>

High Frequency Portable X-ray Unit  
for Large, Small and Mixed Animal Practice



The HF100+ with integrated dual laser pointers is a versatile unit designed to meet the most challenging portable x-ray imaging situations.

Best  
Warranty  
Of Any  
Portable

**Lightest, most powerful 100 kV unit available!**

- **Easily Programmed Memory:**  
*For storage of 5 exposure techniques*
- **mAs/ Time Display Switching Capability**
- **Integrated Dual Laser Pointers:** *Provides for accurate x-ray targeting in bright conditions*
- **High Frequency, Full-wave Rectified Power:**  
*Supplies more radiographic power pound for pound*
- **Automatic Dynamic Line Voltage Compensation:**  
*Unequaled repeatability of results*
- **Exclusive Five-year Limited Warranty:**  
*Worry-free ownership, coverage includes cables*



HF100+ mounted on the 300GS mobile stand for in-clinic imaging.



HF100+ mounted on X100S tripod for field imaging.

MINXRAY<sup>®</sup> INC.  
Leader in High Frequency Portable X-ray Equipment

