

MINXRAY[®] INC.

HF120/60HPPWV **PowerPlus™** **Portable Medical X-ray Unit**



The break-through electronic design of MinXray's new **PowerPlus™** makes it the most powerful portable x-ray unit pound-for-pound available today.

Important features include:

- unique **HighBurst™** electronic design
- up to 120 kVDC
- over 200 mAs
- operation on input of 100-260 V, 50/60 Hz (auto sense & switch)
- CPU control of x-ray tube
- 5 memory stations for storing x-ray exposure techniques
- mAs and time readouts (switchable by toggle)
- low voltage warning light
- compatible with all digital imaging modalities
- seamless integration with existing MinXray stands

MINXRAY[®] INC.
Quality portable x-ray units since 1967

