

The Shield

Patent #
7,082,631

Protect your scanner table from fluid contamination with The Shield. This patented system incorporates a coated foam table pad and reinforced PVC side flaps. The patient rests comfortably on the table pad while the flaps shield the gaps in the table from fluid contamination. The convenience of this system can be adapted to any imaging platform including GE, Siemens, Marconi, Philips, Toshiba and Hitachi.



434-The Shield w/Slide

Benefits of The Shield:

- ScanCoat® Table Pad is fluid resistant and cleanable, creating a sanitary environment, and increasing the life of the pad
- SureChek® Flaps direct fluids away from the patient
- Slide groove technology allows for infinite adjustability of security straps
- Durable Vinyl Security Straps will not tear, while maintaining the patient's position



435-The Shield

Part

Description

434	The Shield w/Slide <i>(Includes Table Pad, Flaps, slide groove runner, Wide Security Strap, Narrow Security Strap)</i>
435	The Shield <i>(Does Not Include Straps or Runner)</i>
436	Replacement Wide Security Strap
437	Replacement Narrow Security Strap



Distributed by: **M.D. McCauley Co., Inc.**
760-C South Rochester Ave.
Ontario, CA 91761 USA

Phone (909) 390-9313
Fax (909) 390-9061

info@xraymdm.com
www.xraymdm.com